



# White Label Products.

Wellness

Skincare

Recovery

Pain relief

Medical

Animal care

and more



**Light Tree**  
TECHNOLOGY.






**Empowering  
Life through  
Technology.**



The global  
medical device  
market is  
expected to  
reach nearly  
\$800 billion  
by 2030.

By 2030 – Global Market Insights

[www.lighttreetechnology.com](http://www.lighttreetechnology.com)



# Table of Contents.

## 1. Introduction to LTT

1.1. About LTT.....	8
1.2. What is White Label? .....	9
1.3. The Process .....	10
1.4. Benefits of working with LTT .....	11

## 2. The Science

2.1. What is Light Therapy?.....	12
2.2. How does it work?.....	12
2.3. Benefits of Light Therapy .....	13

## 3. Beauty & Skin care

3.1. Flexible LED Mask.....	16
3.2. Quad Core LED Silicone Mask.....	17
3.3. LED Silicone Mask .....	18
3.4. LED Hardshell Mask.....	19
3.5. LED Pro Salon Mask .....	20
3.6. Pre Shaped LED Mask 1 .....	21
3.7. Pre Shaped LED Mask 2.....	22
3.8. LED Facial Mask.....	23
3.9. Wired Flat LED Mask 1.....	24
3.10. Wired Flat LED Mask 2 .....	25
3.11. Neck Beauty Care.....	26
3.12. LED Silicone Decollete & Neck Mask .....	27
3.13. LED Silicone Eye Mask.....	28
3.14. Lumina Lift Eye Mask 1.....	29
3.15. Lumina Lift Eye Mask 2.....	30
3.16. Lumina Lift Eye Mask 3 .....	31
3.17. Mermaid Moon Patches.....	32
3.18. Lumera Skin Patches.....	33
3.19. LumiLip System 1.....	34
3.20. LumiLip System 2 .....	35
3.21. LED Eye Mask.....	36



3.22. LED Hair Regrow Cap.....	37
3.23. LED Anti Acne Cleansing Brush.....	38
3.24. LED Beauty Handheld.....	39
3.25. Mini anti-acne handheld.....	40
3.26. Face Contour Applicator .....	41
3.27. Skin Regeneration Guasha .....	42
3.28. Hair Growth Pad .....	43
3.29. LED Beauty Handheld Device.....	44
3.30. LED Anti Acne Pod.....	45
3.31. LED Anti Acne Patches.....	46
3.32. LED Under Eye Patches.....	47
3.33. LED Scalp Massager .....	48
3.34. Anti Acne Patches.....	49

## 4. Recovery & Pain relief

4.1. Universal LED Pad.....	51
4.2. Universal Large LED Pad.....	52
4.3. LED Full Body Pad .....	53
4.4. LED Neck & Upper Back Wrap.....	54
4.5. Shoulder and Upper Back LED Pad with Vibration.....	55
4.6. LED Waist Wrap .....	56
4.7. LED Waist Pad .....	57
4.8. LED Hand Mitt.....	58
4.9. LED Wrist Wrap .....	59
4.10. LED Knee Wrap.....	60
4.11. LED Foot Wrap.....	61
4.12. LED Yoga Mat.....	62
4.13. Large LED Panel .....	63
4.14. Medium LED Panel .....	64
4.15. Small LED Panel .....	65
4.16. LED Slipper .....	66
4.17. LED Power Bed.....	67
4.18. Neck Pillow with Heating .....	68

## 5. Medical


5.1. LED Gum Health Device .....	71
5.2. Oral Care Device .....	72
5.3. LED Cold Sore Wand.....	73
5.4. LED Menstrual Pain Relief Wrap.....	74
5.5. Nasal Accessory.....	75

## 6. Animal care

6.1. LED Relaxation & Pain Relief Pet Bed.....	77
6.2. LED Pet Handmit Brush .....	78
6.3. LED Pet Brush.....	79
6.4. LED Pain Relief Pillow Panel.....	80
6.5. LED Pain Relief Mat.....	81

7. FAQ.....	82
-------------	----

8. Contact.....	86
-----------------	----



LTT is a one-stop shop for Medical and Consumer device development, working with companies of all sizes to co-create innovative products.



# 1. Introduction.

## 1.1. About LTT

Light Tree Technology (LTT) is a global leader in the research, development, and manufacturing of high-quality, technology-driven Medical and Consumer devices.

Our mission is to advance health and wellness through innovative technologies. We accomplish that by researching effective new treatments, designing and developing impactful devices, producing high-quality outputs at competitive prices, and getting regulatory approval for our products. As a one-stop partner for Medical and Consumer device development, LTT collaborates with companies of all sizes to bring innovative products to market.



25

Production Lines

5 million

Annual Output Devices


15.000m<sup>2</sup>

Production Area

## 1.2. What is White Label?

Our white-label services consist of ready-to-sell products we have developed in the past.

We have a full range of products in different categories from Skincare, Wellness, Recovery, Animal care, and Medical. White-label products are a great opportunity to quickly develop, launch, and generate revenue. Our team specializes in quickly transforming a white-label product into your brand's style.



Quickly launch with our ready-to-brand technology-driven devices, offering speed and efficiency to market.



### 1.3. The Process



### 1.4. Benefits of working with LTT

<p><b>15+ years of experience</b> Medical and Consumer devices development and manufacturing.</p>	<p><b>End-to-end R&amp;D process</b> Product concept, design, prototyping, production, and certification.</p>	<p><b>Your innovative partner</b> We create purposeful innovation, working closely with you from idea to execution.</p>
<p><b>Advanced customization</b> Custom solutions for handhelds, masks, wearables and more.</p>	<p><b>25+ New technology concepts per quarter</b> New technology concepts developed and validated each quarter.</p>	<p><b>Cross-industry experience</b> Beauty, wellness, sports, medical, pet care solutions, and more.</p>
<p><b>Certified manufacturing hubs</b> ISO 13485 and FDA/MDSAP-certified facilities in China and India.</p>	<p><b>Partnership-driven execution</b> Partnering for compliance, market fit, and strategic positioning.</p>	




# 2. The Science.

## 2.1. What is Light Therapy?

Light therapy, also known as photobiomodulation, utilizes specific wavelengths of light to stimulate the cellular activity needed for the body's healing processes.

Scientific research has shown that certain light frequencies can promote tissue repair, reduce inflammation, and support circadian rhythm regulation. Its applications range from medical treatments for skin conditions to enhancing cellular activation.



## 2.2. How does it work?

### Cellular Activation

Triggers biochemical reactions in cells similar to the way plants use chlorophyll to convert sunlight into plant tissue.

### Broad-Spectrum Efficacy

Different wavelengths are tailored to target a variety of cellular mechanisms, offering a wide range of therapeutic benefits.

### Biological Regulation

It plays a role in regulating biological functions, which may include hormonal balance and immune system response.

## 2.3. Benefits of Light Therapy

Light is just like food and water. Our body needs it to survive and be healthy.

Light Therapy was discovered and introduced by NASA in the mid 20th century. Initially, it was primarily used for astronauts traveling to space to ensure their health and to address certain physical and mental issues.

The cornerstone of why light therapy works is because it increases blood circulation and ATP production, thereby providing more energy for your body's cells.

These two benefits form the foundation for the majority of treatments. Additional benefits are mainly treatment-specific.

### Treatments

- Anti-Acne
- Anti-Aging
- Oral Care
- Hair Growth
- Pain Relief
- Muscle Recovery
- Pet Care
- Disinfection
- Teeth Whitening
- Other Treatments



Elevate your wellness brand with LTT's versatile LED therapy products. Expand into skin care, hair growth, oral care, and muscle recovery with ease.

Join the wellness movement and diversify your offerings today.



**Skin care**




**Oral care**




**Hair care**



**Muscle recovery**



## 3. Beauty & Skin care.



LED light therapy is clinically-proven and offers a therapeutic treatment that is safe, non-invasive, and chemical-free. LED-masks have been shown to be effective in clearing acne, reducing wrinkles, evening skin tones, reducing redness, and revitalizing, soothing, and hydrating skin.

### 3.1. Flexible LED Mask

We have modernized the traditional concept of a face mask by using current light therapy technology to create something special.

The LED Mask combines benefits from various light frequencies, including infrared and its skin boosting abilities. This combination penetrates deep into the skin and offers the most diverse solution for skin imperfections.



#### Benefits:

- **Skin boosting**  
Two powerful at home treatments.
- **100% UV Free**  
No harmful UV rays and no recovery time.
- **Non-invasive**  
Safe for all type of skin.
- **Chemical free**  
No need for abrasive chemical products.
- **Dutch design**  
We strive for quality and perfection. All products meet high quality standards and are tested before leaving our factory.
- **Clinically tested and proven**  
Our products are clinically tested and proven to be safe and effective.



#### One solution for all skin issues

- Wrinkles
- Eyes bags
- Dark circles
- Puffiness
- Large pores
- Dull skin
- Sensitive acne

### 3.2. Quad Core LED Silicone Mask

Product code: MK-03

#### Size of product

L706 x W227 x H5 (mm)



#### Weight of product

370g

#### Specifications

<b>Amount of LEDs</b>	66 triple-core blue + red + infrared LEDs
<b>Wavelength(s)</b>	465nm (blue) 620nm (red) 850nm (infrared)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Adjustable strap, manual, controller, USB-C charging cable, eye shield





### 3.3. LED Silicone Mask


Product code: MK-66-RGB

**Size of product**  
L310x W223 x H5 (mm)

**Weight of product**  
152g

#### Specifications

<b>Amount of LEDs</b>	66 triple-core red + blue LEDs 20 single-core infrared LEDs
<b>Wavelength(s)</b>	465nm (blue) 620nm (red) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Safety goggles, power adapter, adjustable strap, manual, USB-C charging cable



### 3.4. LED Hardshell Mask

Product code: MJ-09

**Size of product**  
L166 x W198 x H108 (mm)

**Weight of product**  
800g

#### Specifications

<b>Amount of LEDs</b>	50 dual-core red + infrared LEDs 30 dual-core red + blue LEDs
<b>Wavelength(s)</b>	415nm (blue) 633nm (red) 830nm (infrared)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup> 34 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C cable, head strap, built-in eye protection, manual, luxury pouch



### 3.5. LED Pro Salon Mask

Product code: MK-177

#### Size of product

H180 (cm)

#### Weight of product

1050g

#### Specifications

##### Amount of LEDs

Eyes	30 triple-core red + blue + green LEDs
Nose	12 triple-core red+ blue + green LEDs
Face & forehead	135 triple-core red+ blue + green LEDs
Lips	20 single-core blue LEDs
Full face	70 single-core infrared LEDs

##### Wavelength(s)

400nm (purple)  
440nm (blue)  
500nm (cyan)  
520nm (green)  
565nm (yellow)  
590nm (orange)  
625nm (red)

##### Treatment time

No standard time

##### Frequency of use

Daily

##### Type of treatment

Beauty & skin care

##### Power




30 mW/cm<sup>2</sup>

##### Charging method

Direct wall socket

##### Accessories

Charging cable



### 3.6. Pre Shaped LED Mask 1

Product code: MJ-66B

#### Specifications

Amount of LEDs	66 pcs triple LEDs
Wavelength(s)	415nm 630nm 830nm
Treatment time	10 minutes
Type of treatment	Beauty & skin care
Power	30 mW/cm <sup>2</sup>
Accessories	Straps, USB cable, manual, package

### 3.7. Pre Shaped LED Mask 2

Product code: MK-66F

#### Specifications

<b>Amount of LEDs</b>	78 pcs triple LEDs
<b>Wavelength(s)</b>	415nm 630nm 830nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, USB cable, manual, package





### 3.8. LED Facial Mask

Product code: MK-78E

#### Specifications

<b>Amount of LEDs</b>	78 pcs triple LEDs
<b>Wavelength(s)</b>	415nm 630nm 830nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, controller, USB cable, manual, package




### 3.9. Wired Flat LED Mask 1

Product code: MK-180

#### Specifications

<b>Amount of LEDs</b>	78 pcs triple LEDs
<b>Wavelength(s)</b>	415nm 630nm 830nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, Controller, USB cable, manual, package




### 3.10. Wired Flat LED Mask 2

Product code: MK-07

#### Specifications

<b>Amount of LEDs</b>	78 pcs triple LEDs
<b>Wavelength(s)</b>	415nm 630nm 830nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, Controller, USB cable, manual, package




### 3.11. Neck Beauty Care

Product code: ZLD-80

#### Specifications

<b>Amount of LEDs</b>	49 pcs triple LEDs
<b>Wavelength(s)</b>	630nm 830nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, controller, USB cable, manual, package



### 3.12. LED Silicone Decollete & Neck Mask

Product code: ZLD-50C

#### Size of product

L340 x W210 x H5 (mm)

#### Weight of product

268g


#### Specifications

<b>Amount of LEDs</b>	49 dual-core red + infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 830nm (infrared)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Long and short strap, controller, USB-C charging cable, manual



### 3.13. LED Silicone Eye Mask

Product code: EY-20



**Size of product**

L299 x W77 x H5 (mm)

**Weight of product**

92g

**Specifications**

<b>Amount of LEDs</b>	20 dual-core red + amber LEDs 20 dual-core red + infrared LEDs
<b>Wavelength(s)</b>	605nm (amber) 630nm (red) 660nm (red) 880nm (infrared)
<b>Treatment time</b>	3 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Adjustable strap, manual, USB-C charging cable




### 3.14. Lumina Lift Eye Mask 1

Product code: EY-40R

**Specifications**

<b>Amount of LEDs</b>	40 pcs double LEDs
<b>Wavelength(s)</b>	605nm 633nm 660nm 830nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, USB cable, manual, package





### 3.15. Lumina Lift Eye Mask 2

Product code: EY-20DA

#### Specifications

<b>Amount of LEDs</b>	40 pcs double LEDs
<b>Wavelength(s)</b>	605nm 633nm 660nm 830nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, USB cable, manual, package




### 3.16. Lumina Lift Eye Mask 3

Product code: EY-40F

#### Specifications

<b>Amount of LEDs</b>	40 pcs double LEDs
<b>Wavelength(s)</b>	605nm 633nm 660nm 830nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65 mW/cm <sup>2</sup>
<b>Accessories</b>	Straps, USB cable, manual, package





### 3.17. Mermaid Moon Patches

Product code: MT-47

#### Specifications

<b>Amount of LEDs</b>	6 pcs double LEDs
<b>Wavelength(s)</b>	460nm 630nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	5~10mW/cm <sup>2</sup>
<b>Accessories</b>	Facial sticker, USB cable, manual, package






### 3.18. Lumera Skin Patches

Product code: MT-43

#### Specifications

<b>Amount of LEDs</b>	6 pcs double LEDs
<b>Wavelength(s)</b>	415nm 630nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	Mini controller, facial sticker, USB cable, manual, package





### 3.19. LumiLip System 1

Product code: OR-38

#### Specifications

<b>Amount of LEDs</b>	38 pcs double LEDs
<b>Wavelength(s)</b>	605nm 630nm 660nm 880nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65W/cm <sup>2</sup>
<b>Accessories</b>	USB cable, manual, package




### 3.20. LumiLip System 2

Product code: OR-01

#### Specifications

<b>Amount of LEDs</b>	38 pcs double LEDs
<b>Wavelength(s)</b>	605nm 630nm 660nm 880nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	65 mW/cm <sup>2</sup>
<b>Accessories</b>	Charging case, USB cable, manual, package





### 3.21. LED Eye Mask

Product code: T26-EY

#### Size of product

L310 x W80 x H5 (mm)

#### Weight of product

56g

#### Specifications

<b>Amount of LEDs</b>	26 dual-core red and infrared LEDs
<b>Wavelength(s)</b>	625nm (red) 850nm (infrared)
<b>Treatment time</b>	3 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	20 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual, strap



### 3.22. v

Product code: T-120-HRC

#### Size of product

L160 x W160 x H200 (mm)

#### Weight of product

300g

#### Specifications

<b>Amount of LEDs</b>	120 single-core red LEDs
<b>Wavelength(s)</b>	640nm (red)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	20 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Manual, USB-C charging cable



### 3.23. LED Anti Acne Cleansing Brush

Product code: HD-05-I

#### Size of product

L155 x W50 x H25 (mm)

#### Weight of product

48g

#### Specifications

<b>Amount of LEDs</b>	5 dual-core blue + red LEDs
<b>Wavelength(s)</b>	415nm (blue) 660nm (red)
<b>Treatment time</b>	3 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Manual, USB-C charging cable



### 3.24. LED Beauty Handheld

Product code: SB-01N

#### Size of product

L137 x W61 x H28 (mm)



#### Weight of product

345g

#### Specifications

<b>Amount of LEDs</b>	9 single-core red LEDs 9 single-core green LEDs 9 dual-core infrared & yellow LEDs 9 dual-core infrared & blue LEDs
<b>Wavelength(s)</b>	470nm (blue) 520nm (green) 590nm (yellow) 650nm (red) 940nm (infrared)
<b>Treatment time</b>	5, 10 or 15 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup> 35 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual





### 3.25. Mini anti-acne handheld

Product code: M-03

#### Specifications

<b>Amount of LEDs</b>	3 pcs double LEDs
<b>Wavelength(s)</b>	415nm
<b>Treatment time</b>	3 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	20 - 65W/cm <sup>2</sup>
<b>Accessories</b>	USB cable, manual, package




### 3.26. Face Contour Applicator

Product code: HD-44

#### Specifications

<b>Amount of LEDs</b>	3 pcs double LEDs
<b>Wavelength(s)</b>	630±10nm
<b>Treatment time</b>	5 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30 mW/cm <sup>2</sup>
<b>Accessories</b>	USB cable, manual, package





### 3.27. Skin Regeneration Guasha

Product code: HD-186

#### Specifications

<b>Amount of LEDs</b>	3 pcs double LEDs
<b>Wavelength(s)</b>	630±10nm
<b>Treatment time</b>	5 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	30W/cm <sup>2</sup>
<b>Accessories</b>	USB cable, manual, package



### 3.28. Hair Growth Pad

Product code: MP12

#### Specifications

<b>Amount of LEDs</b>	130 pcs Double LEDs:
<b>Wave</b>	650nm
<b>Treatment time</b>	10 minutes
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	5 mW/cm <sup>2</sup>
<b>Accessories</b>	Controller, USB cable, manual, package





### 3.29. LED Beauty Handheld Device

Product code: HD-03AL

**Size of product**

L173 x W74 x H15 (mm)

**Weight of product**

133g

**Specifications**

<b>Amount of LEDs</b>	12 triple-core red + infrared + violet LEDs
<b>Wavelength(s)</b>	415nm (violet) 630nm (red) 880nm (infrared)
<b>Treatment time</b>	3 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	60 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



### 3.30. LED Anti Acne Pod

Product code: MT-15FG

**Size of product**

L72 x W44 x H30 (mm)

**Weight of product**

24g

**Specifications**

<b>Amount of LEDs</b>	15 single-core violet LEDs
<b>Wavelength(s)</b>	415nm (violet)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Beauty & skin care
<b>Power</b>	5 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, 60 strips, manual



### 3.31. LED Anti Acne Patches

Product code: MT-10C

#### Size of product

L86 x W41 x H5 (mm)

#### Weight of product

10g per patch

#### Specifications

**Amount of LEDs** 12 single-core blue LEDs

12 single-core red LEDs

**Wavelength(s)** 415nm (blue)

630nm (red)

**Treatment time** 3 minutes

**Frequency of use** Daily

**Type of treatment** Beauty & skin care

**Power** 30 mW/cm<sup>2</sup>

**Charging method** USB-C

**Accessories** 2 patches, USB-C charging cable with integrated controller, 1 travel pouch, 120 strips, manual



### 3.32. LED Under Eye Patches

Product code: MT-10D

#### Size of product

L86 x W41 x H5 (mm)

#### Weight of product

10g per patch

#### Specifications

**Amount of LEDs** 11 dual-core red & blue LEDs (per patch)

**Wavelength(s)** 415nm (blue)

630nm (red)

**Treatment time** 3 minutes


**Frequency of use** Daily

**Type of treatment** Beauty & skin care

**Power** 30 mW/cm<sup>2</sup>

**Charging method** USB-C

**Accessories** 2 beauty patches, USB-C charging cable with integrated controller, 1 travel pouch, 120 strips, manual



### 3.33. LED Scalp Massager

Product code: SM-12

#### Size of product

L90 x W80 x H80 (mm)

#### Weight of product

73g

#### Specifications

**Amount of LEDs** 12 single-core red LEDs

**Wavelength(s)** 660nm (red)

**Treatment time** 3 minutes


**Frequency of use** Daily

**Type of treatment** Beauty & Hair care

**Power** 20 mW/cm<sup>2</sup>

**Charging method** USB-C

**Accessories** USB-C charging cable,  
manual



### 3.34. Anti Acne Patches

Product code: MT-10C



#### Size of product

L51 x W36 x 14(mm)

#### Weight of product

14g

#### Specifications

**Amount of LEDs** 6 LEDs

**Wavelength(s)** 630nm (red)  
415nm (blue)

**Treatment time** 3 minutes

**Frequency of use** Daily

**Type of treatment** Anti acne

**Power** 9.5 mW/cm<sup>2</sup>  
21 mW/cm<sup>2</sup>

**Charging method** USB-C

**Accessories** Wireless with charging dock





# 4. Recovery & Pain relief.

Clinically proven to enhance performance and accelerate recovery, LED light therapy has become an indispensable tool for athletes, both professional and non-professional, seeking to optimize their physical well-being. From relieving pain and facilitating faster recovery post-exercise to improving circulation and promoting overall vitality, LED light therapy offers a scientifically backed approach to sports recovery and rejuvenation.

## 4.1. Universal LED Pad

Product code: ZLD-60



### Size of product

L270 x W140 x H5 (mm)

### Weight of product

120g

### Specifications

<b>Amount of LEDs</b>	30 single-core Red LEDs 30 single-core Infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	25 mW/cm <sup>2</sup> 60 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, power adapter, user manual, long & short strap



## 4.2. Universal Large LED Pad

Product code: ZLD-160



### Size of product

L380 x W190 x H5 (mm)

### Weight of product

250g


### Specifications

<b>Amount of LEDs</b>	80 single-core red LEDs 80 single-core infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 660nm (red) 810nm (infrared) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	25 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



## 4.3. LED Full Body Pad

Product code: ZLD-320



### Size of product

L272.8 x W457.2 x H5 (mm)

### Weight of product

500g


### Specifications

<b>Amount of LEDs</b>	160 single-core red LEDs 160 single-core infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 660nm (red) 810nm (infrared) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	10-20 mW/cm <sup>2</sup> 30-50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



#### 4.4. LED Neck & Upper Back Wrap

Product code: ZLD-60-NW



##### Size of product

L700 x W49 x 5 (mm)

##### Weight of product

460g

##### Specifications

<b>Amount of LEDs</b>	30 single-core red LEDs 30 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	25 mW/cm <sup>2</sup> 60 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, power adapter, user manual



#### 4.5. Shoulder and Upper Back LED Pad with Vibration

Product code: ZLD-60-BO



##### Size of product

L700 x W236 x 5 (mm)

##### Weight of product

350g

##### Specifications

<b>Amount of LEDs</b>	30 single-core red LEDs 30 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	40 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, power adapter, user manual



#### 4.6. LED Waist Wrap

Product code: ZLD-60-BP



##### Size of product

L1200 x W240 x H5 (mm)

##### Weight of product

340g

##### Specifications

<b>Amount of LEDs</b>	30 single-core red LEDs 30 single-core Infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	25 mW/cm <sup>2</sup> 60 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, power adapter, user manual



#### 4.7. LED Waist Pad

Product code: HY-60Q



##### Size of product

L745 x W497 x 22.5 (mm)

##### Weight of product

800g

##### Specifications

<b>Amount of LEDs</b>	30 single-core red LEDs 30 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	7 days /week for 6 weeks
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	40±10 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	Manual, USB-C charging cable



#### 4.8. LED Hand Mitt

Product code: ST-40

##### Size of product




L260 x W150 (mm)

##### Weight of product

138g

##### Specifications

<b>Amount of LEDs</b>	20 single-core red LEDs 20 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	30 mW/cm <sup>2</sup> 35 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



#### 4.9. LED Wrist Wrap

Product code: HW-24

##### Size of product

L420 x W320 x H80 (mm)

##### Weight of product

160g

##### Specifications

<b>Amount of LEDs</b>	12 single-core red LEDs 12 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	50mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual

#### 4.10. LED Knee Wrap

Product code: HX-96



##### Size of product

L620 x W295 x H5 (mm)

##### Weight of product

205g


##### Specifications

<b>Amount of LEDs</b>	48 single-core red LEDs 48 single-core infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 660nm (red) 810nm (infrared) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	20 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



#### 4.11. LED Foot Wrap

Product code: HJ-38H



##### Size of product

L575 x W410 x H5 (mm)

##### Weight of product

150g


##### Specifications

<b>Amount of LEDs</b>	20 single-core red LEDs 18 single-core infrared LEDs
<b>Wavelength(s)</b>	650nm (red) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



#### 4.12. LED Yoga Mat

Product code: BL-400-YM



##### Size of product

L820 x W610 x H8 (mm)

##### Weight of product

1700g


##### Specifications

<b>Amount of LEDs</b>	400 dual-core red & infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 850nm (infrared)
<b>Treatment time</b>	No standard time
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	20 mW/cm <sup>2</sup> 25 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct to wall socket
<b>Accessories</b>	Manual, charging cable, remote with light intensity and timer



#### 4.13. Large LED Panel

Product code: EA-150



##### Size of product

L945 x W221 x H59.5 (mm)

##### Weight of product

6300g

##### Specifications

<b>Amount of LEDs</b>	30 dual-core red LEDs 120 dual-core infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 650nm (red) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily use
<b>Type of treatment</b>	Care
<b>Power</b>	180-210 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct wall socket
<b>Accessories</b>	Manual, charging cable, goggles



#### 4.14. Medium LED Panel

Product code: EA-60

##### Size of product


L387 x W221 x H59.5 (mm)

##### Weight of product

3500g

##### Specifications

<b>Amount of LEDs</b>	10 dual-core red LEDs 50 dual-core infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 650nm (red) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	180-210 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct to wall socket
<b>Accessories</b>	Manual, charging cable, goggles



#### 4.15. Small LED Panel

Product code: EA-12

##### Size of product


L96 x W137 x H38 (mm)

##### Weight of product

1500g

##### Specifications

<b>Amount of LEDs</b>	12 dual-core red & infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Recovery & pain relief
<b>Power</b>	100 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct to wall socket
<b>Accessories</b>	Manual, charging cable, goggles



#### 4.16. LED Slipper

Product code: XZ-60

##### Size of product

L324.3 x W133 (mm)

##### Weight of product

200g

##### Specifications

**Amount of LEDs** 60 single core LED

**Wavelength(s)** 660nm (red)  
880nm (infrared)

**Treatment time** 20 minutes

**Frequency of use** Daily

**Type of treatment** Recovery  
& pain relief

**Power** 50 +/-5 mW/cm<sup>2</sup>  
50 +/-5 mW/cm<sup>2</sup>

**Charging method** USB-C

**Accessories** Charging cable,  
manual, adapter



#### 4.17. LED Power Bed

Product code: EAS24-P

##### Size of product

L189 x W77 x H118 (cm)

##### Weight of product

166kg

##### Specifications

**Amount of LEDs** 3840 quad-core LEDs

**Wavelength(s)** 485nm (blue)  
660nm (red)  
760nm (red)  
880nm (infrared)

**Treatment time** No standard time

**Frequency of use** Daily

**Type of treatment** Recovery  
& pain relief

**Power** 90-110m W/cm<sup>2</sup>

**Charging method** Direct to wall socket

#### 4.18. Neck Pillow with Heating

Product code: ZT-24

##### Size of product

L285 x W260 x H140 (mm)

##### Weight of product

420g

##### Specifications

Amount of LEDs	24 dual-core red + infrared LEDs
Wavelength(s)	660nm (red) 880nm (infrared)
Treatment time	20 minutes
Frequency of use	Daily
Type of treatment	Recovery & pain relief
Power	27 +/-5 mW/cm <sup>2</sup> 42 +/-5 mW/cm <sup>2</sup>
Charging method	USB-C
Accessories	Charge cable, manual, adapter



FDA



# 5. Medical.



Backed by extensive research, LED light therapy has been shown to reduce bone, muscle, and joint pain, providing relief to individuals suffering from conditions such as arthritis, temporomandibular joint disorder (TMJD), and sports injuries. Additionally, LED light therapy has recently emerged as a game-changer in dental care, offering effective solutions for oral health issues such as gum disease, inflammation, and oral infections.

## 5.1. LED Gum Health Device

Product code: YT-08

### Size of product

L85 x W65 x H25 (mm)

### Weight of product

37g

### Specifications

<b>Amount of LEDs</b>	16 dual-core red + infrared LEDs
<b>Wavelength(s)</b>	630nm (red) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Medical
<b>Power</b>	20 mW/cm <sup>2</sup>
<b>Charging method</b>	Magnetic connector to USB-A
<b>Accessories</b>	Charging cable, manual, product case



## 5.2. Oral Care Device

Product code: YT-11

### Specifications

Amount of LEDs	20 pcs double LEDs
Wavelength(s)	460nm 650nm
Treatment time	8 minutes
Type of treatment	Medical
Power	20 mW/cm <sup>2</sup>
Accessories	Magnetic cable, manual, package



## 5.3. LED Cold Sore Wand

Product code: HD-27

### Specifications

Amount of LEDs	2 single-core red & blue LEDs
Wavelength(s)	415nm (blue) 630nm (red)
Treatment time	3 minutes
Frequency of use	Daily
Type of treatment	Medical
Power	30 mW/cm <sup>2</sup>
Charging method	USB-C
Accessories	USB-C charging cable, manual



## 5.4. LED Menstrual Pain Relief Wrap

Product code: ZLD-80-MP

### Size of product


L320 x W170 x H5 (mm)

### Weight of product

130g

### Specifications

<b>Amount of LEDs</b>	80 dual-core red and infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 880nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Medical
<b>Power</b>	20 mW/cm <sup>2</sup> 35 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual, strap



## 5.5. Nasal Accessory

Product code: BZ-02

### Size of product


L42.7 x W13.7 x H28 (mm)

### Weight of product

20g

### Specifications

<b>Amount of LEDs</b>	2 dual-core red & blue LEDs
<b>Wavelength(s)</b>	415nm (blue) 660nm (red)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Once a day, 3-5 times a week
<b>Type of treatment</b>	Medical
<b>Power</b>	35 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct to PD 12V DC charger
<b>Accessories</b>	Type-C to Type-C cable, silicone caps, metal box, color box, user manual






# 6. Animal care.



When your animal suffers from joint pain, wounds, or needs post-surgery care, light therapy can offer relief. In animal rehabilitation, red light therapy is a commonly proven method for various conditions. It helps reduce pain and inflammation and aids in wound care and healing. With its many benefits, red light therapy has become a go-to treatment for many canine and equine conditions.

## 6.1. LED Relaxation & Pain Relief Pet Bed



### Specifications

<b>Amount of LEDs</b>	30 single-core red LEDs 30 single-core infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 810nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Animal
<b>Power</b>	25 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C or direct to wall socket
<b>Accessories</b>	Charging cable, manual



## 6.2. LED Hand Mitt

### Specifications

<b>Amount of LEDs</b>	22 dual-core red and infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 810nm (infrared)
<b>Treatment time</b>	5 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Animal
<b>Power</b>	25 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C or direct to wall socket
<b>Accessories</b>	Charging cable, manual



## 6.3. LED Pet Brush

### Specifications

<b>Amount of LEDs</b>	4 dual-core red and infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 810nm (infrared)
<b>Treatment time</b>	3 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Animal
<b>Power</b>	25 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C
<b>Accessories</b>	USB-C charging cable, manual



## 6.4. LED Pain Relief Pillow Panel

### Specifications

<b>Amount of LEDs</b>	30 single-core red and infrared LEDs
<b>Wavelength(s)</b>	660nm (red) 850nm (infrared)
<b>Treatment time</b>	10 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Animal
<b>Power</b>	25 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	USB-C or direct to wall socket
<b>Accessories</b>	Charging cable, manual




## 6.5. LED Pain Relief Mat

### Specifications

<b>Amount of LEDs</b>	320 single-core red LEDs
<b>Wavelength(s)</b>	630nm (red) 660nm (red) 810nm (infrared) 830nm (infrared) 850nm (infrared)
<b>Treatment time</b>	20 minutes
<b>Frequency of use</b>	Daily
<b>Type of treatment</b>	Animal
<b>Power</b>	25 mW/cm <sup>2</sup> 50 mW/cm <sup>2</sup>
<b>Charging method</b>	Direct to wall socket
<b>Accessories</b>	Charging cable, manual

# 7. FAQ.



## LIGHT THERAPY

### 1. What is Light Therapy?

Light therapy, also known as phototherapy, is a (medical) treatment that involves exposure to specific wavelengths of light to treat various conditions.

### 2. How does Light Therapy work?

Light therapy works by stimulating the mitochondria, or the powerhouse, of your cells. This cellular stimulation can result in an increased production of cellular energy. With this extra energy, your cells may be able to regenerate and repair themselves faster.

### 3. What conditions can Light Therapy treat?

It's effective for various conditions including skin issues like acne and signs of aging, mood disorders like SAD, and pain relief for muscle and joint discomfort.

### 4. Is Light Therapy safe?

Yes, when used as directed, light therapy is safe. However, it's important to follow guidelines and consult with a healthcare professional, especially for individuals with certain health conditions.

### 5. How long does a Light Therapy session last?

The duration of a light therapy session can vary depending on the condition being treated and the specific device used. For example, an anti-aging treatment can take from three to twenty minutes depending on what kind of device you are using.

### 6. Can Light Therapy be used at home?

Yes, there are home-use light therapy devices available like a handheld or facial LED mask.

### 7. Are there any side effects of Light Therapy?

Side effects are typically minimal but can include temporary discomfort, such as eyestrain, headache, or nausea.

### 8. Is Light Therapy clinically proven?

Yes, there are many clinical studies proving the different elements of light therapy.

### 9. Is Light Therapy FDA cleared or approved?

Yes, several treatments are FDA cleared including:

- Oral care
- Anti-acne
- Anti-aging
- Pain relief
- Muscle recovery

### 1. What is Light Tree Technology's primary focus?

Light Tree Technology focuses on the research, development, and manufacturing of technology-driven Medical and Consumer devices. Our core expertise includes LED light therapy and other evidence-based technologies used across skincare, wellness, medical, and animal care applications.

### 2. Where are Light Tree Technology's products and services available?

We operate across the globe. Headquartered in the Netherlands, we also have factories in China and India. We have a secondary office in the USA which also contains a fulfillment location.

### 3. What makes Light Tree Technology unique?

Light Tree Technology is a one-stop technology partner for brands developing medical and consumer devices. Whether you want to go for Private or White label, we can help you without the involvement of other suppliers or companies. At the moment, we are working for clients in the skincare, wellness, medical, and animal care industries.

### 4. How does Light Tree Technology ensure the quality and efficacy of its products?

All Light Tree Technology locations are ISO13485 and MDSAP certified which means we have to follow certain processes that lead to higher quality.

### 5. Can individuals and businesses collaborate or partner with Light Tree Technology?

Yes, we have helped both businesses and individuals create a product to bring to market. We are here to support you.

### 6. Who is Light Tree Technology?

Light Tree Technology specializes in Research, Development and Manufacturing of technology-driven Medical and Consumer devices. We are a one-stop shop for Private and White Label services with offices in the Netherlands, China, India, and the USA.

### 7. What services does Light Tree Technology offer?

Private label and White label.

### 8. Which industries do you operate in?

We support clients within the beauty, medical, health, sports, and wellness industries.

### 1. What range of services does Light Tree Technology offer?

Light Tree Technology mainly offers White and Private label services. Besides that our company offers additional services such as Intellectual Property, Fulfillment and Marketing support.

### 2. Can Light Tree Technology customize Light Therapy treatments for individual needs?

Yes, this is possible. Our team of professionals is here to help guide your vision and requirements into a final product.

### 3. What is the difference between Private and White label?

The main difference between Private and White label is "Customisation". If you choose Private label, we will develop your product from scratch. When you choose White label, you can choose from a catalogue or ready to sell products with less customization.

### 4. How much is the MOQ for White label?

Typically 500.

### 5. How do you involve your clients in the design and development process?

First of all, the client determines the product requirements. The product will be designed and developed according to those requirements. During the design and development process, the client has several opportunities to confirm the design. We will produce a prototype that will also require your final approval.

### 6. What is White label?

Selling a device that is ready-to-go with your own logo and colours. Packaging can also be changed.

# 8. Contact.

Ready to  
develop your own  
light therapy device?

Feel free  
to contact us!

**Call Us**

WhatsApp: +17866471640

**Email**

hello@lighttreeventures.com

**HQ - The Netherlands**

Light Tree Ventures Europe B.V.  
Saturnusstraat 95  
2516 AG, The Hague  
The Netherlands

**USA**

Light Tree Ventures USA LLC  
6240 NE 4th CT  
33138, Miami, FL  
The United States of America

**China**

Shenzhen Light Tree Ventures  
Medical Technology Co. Ltd.  
Building 2, Huaide Cuihu Industrial  
Park, Xintang Road, Bao'an District  
518103, Shenzhen  
China

**India**

Light Tree Technology India Pvt. Ltd.  
Plot No. GP-47, Sector 18, Maruti  
Udyog  
122008, Gurugram, Haryana  
India





Light Tree Technology Europe BV  
Saturnusstraat 95, Office 2.01  
2516 AG, The Hague  
The Netherlands  
[www.lighttreetechnology.com](http://www.lighttreetechnology.com)